

ACTIN-TOOLKITS

Actin-Based Bioassays for
Functional & Structural Ligand Analysis



Analytical Biochemistry

Molecular Cell Biology

Proteomics

Structural Biology

Molecular Medicine

...



www.hypermol.com

Contents

Kit Contents	2
Shipping, Storage and Handling Conditions	2
Product Warranty	2
Preface	3
Introduction to Actin and the Actin Cytoskeleton	4
<i>In Vitro</i> Polymerization of Actin	5
Experimental Conditions	6
Protocol Section	
Protocol 1. Preparation G-Actin	7
Storage of G-Actin	7
Protocol 2. Preparation F-Actin	8
Storage of F-Actin	8
Protocol 3. F-Actin Cosedimentation at High Speed	9
Protocol 4. F-Actin Bundling Assay at Low Speed	11
Troubleshooting Guide	13
Ordering Information	15
Hypermol Services	16

Notice to Purchaser:

For Use in Research Only (RUO).

For *In Vitro* Use Only.

Not for Use in Human or Animal Diagnostic or Therapeutic Processes.

Handbook Actin-Toolkit: F-Actin Binding

© by HYPERMOL EK, all rights reserved.

All products from HYPERMOL are made in Germany.

Kit Content

Actin-Toolkit F-Actin Binding

Cat.# 8010-01 / Cat.# 8100-01

Actin	4x1mg
PolyMix ¹ , PolyMix6 ² , PolyMix8 ³	1x1ml (10x stock sol.) each
MonoMix ⁴	4x50ml
MgCl ₂ (1M)	1x500μl
Handbook	1x

¹ 1M KCl, 0.02M MgCl₂, 0.01M ATP, 0.1M imidazole pH 7.4

³ 1M KCl, 0.02M MgCl₂, 0.01M ATP, 0.1M imidazole pH 8.0

² 1M KCl, 0.02M MgCl₂, 0.01M ATP, 0.1M MES pH 6.0

⁴ 0.1mM CaCl₂, 0.5mM DTT, 0.4mM ATP, 2mM Tris pH 8.2

Shipping, Storage and Handling Conditions

Actin-Toolkits are shipped at ambient temperature. Proteins and ATP-containing buffers are stable for at least 3 months when frozen at -20°C , and for at least 6 months when frozen at -70°C upon arrival. Avoid repeated freeze/thaw of ATP-containing components and of proteins. Solubilized proteins and ATP-containing solutions must be kept on ice and used as described in the Protocol Section.

Product Warranty

Hypermol guarantees the quality and product performance described in this handbook only when products are frozen upon arrival as mentioned above. We do not take any guarantee for uses of our products other than described here. This product is designed for *in vitro* use only and not for use in diagnostic processes.

Should any product fail to perform as guaranteed due to reasons other than misuse or should not meet your expectations, please first contact our TechnicalService (techserv@hypermol.com or ++495219876230) within 5 working days and then return the product to Hypermol as advised. We reserve the right to test the performance of returned products in order to suggest replacement free of charge or refund of the purchase price. The buyer obtains a copy of our "Terms and Conditions of Sale" before ordering, and agrees to this by ordering.

Preface

Identification of an actin-binding protein is a major step forward in research. Today about two hundred proteins are known to either possess direct or cryptic binding sites for actin.

- Actin-Toolkits were developed to safely guide experiments with actin of highest quality.
- The handbook provides background information and protocols for successful and error-free handling.
- Actin-Toolkit proteins are fully biologically active.
- Actin-Toolkits are user-friendly all-in-one applications.

Actin Toolkits are invaluable tools to analyze the biological activity of protein ligands, especially for recombinant proteins, fragments or mutants. Identification and mapping of actin-binding sites in full length proteins or fragments are examples for the use of these assays.

Introduction to Actin and the Actin Cytoskeleton

Actin is one of the most abundant proteins of eukaryotic cells. Comprising 5 to 10% of the total cellular protein, actin turned out to be a key protein of cellular architecture and thus keeper of cellular functions.

The most characteristic feature of single actin molecules (globular or G-actin, Mr=42kD, single polypeptide chain) is to polymerize into double helical filaments (filamentous or F-actin) of several micrometres length *in vitro* and *in situ*. This polymerization process is readily initiated at physiological salt concentrations in the presence of ATP.

Today more than two hundred proteins are known to possess one or more actin binding sites. Some proteins can readily bind to actin; some have to undergo ligand induced conformational changes to bind to actin. Several dozens of proteins directly modulate either the state or the conformation of F- or G-actin.

In addition to the filamentous actin incorporated into the cytoskeleton, cells have a rather variable pool of unpolymerized actin (30-50% of the total actin). Actin sequestering proteins like thymosin β 4 take control of the G-actin pool, which would otherwise polymerize considering intracellular conditions.


Apart from these direct actions on actin - like polymerization, nucleation, capping, depolymerization, severing, bundling etc. - the function of many actin-binding proteins is to support the different states of the actin-cytoskeleton in order to follow the demands of the cellular life (movement, cell division, signalling, etc.).

Due to its central function for the cell, it is not surprising, that the amino acid sequence of actin is highly conserved throughout evolution. Human skeletal muscle actin is identical in sequence to muscle actin in mouse, rat, rabbit, chicken, and beef.

In Vitro Polymerization of Actin

An increase of the salt concentration of a low salt buffer (e.g. <math><5\text{mM}</math> KCl) in the presence of 0.1mM ATP to a nearly intracellular salt concentration ($\sim 100\text{mM}$ KCl, 2mM MgCl_2) induces the polymerization of G-actin to F-actin. In biochemical terms, the polymerization of actin is a non-covalent association of actin monomers. Each actin monomer (G-actin) has two actin binding sites and is incorporated into the double helical filament via head-to-tail association. The result of this polymerization is a polar, right handed double-stranded actin helix (synonym. actin filament) with a pitch of 73.6nm per turn and a diameter of $\sim 8\text{nm}$.

The polymerization of actin is divided into three steps. During the first step (nucleation), actin monomers form an unstable dimer, which becomes stabilized by the addition of another G-actin molecule to form a trimer. Actin-trimers strongly favour the elongation of the actin filament by further monomer addition. Both, temperature and concentration alter polymerization kinetics.



G-actin is an ATPase converting ATP in $\text{ADP} \cdot \text{P}_i$. Either ATP or ADP is bound in the nucleotide binding pocket of the actin molecule. As the affinity of actin is higher for ATP than for ADP, the ADP is exchanged. Free actin monomers hydrolyse ATP

at a very low rate. Upon polymerization the process of hydrolysis is significantly accelerated.

In the early stage of elongation the polarity of actin filaments becomes obvious. Actin polymerizes from both ends of the filament, but the rate of polymerization at either end is different. A fast growing end (plus-end or barbed end) and a slowly growing minus-end (pointed end) are distinguished. The polarity of actin filaments is the consequence of the head-to-tail association of the monomers. The critical concentration (C_c) for polymerization is $\sim 0.1\mu\text{M}$ for the plus-end and $\sim 0.6\mu\text{M}$ for the minus-end. As a consequence, addition of monomers to the plus-end occurs below the C_c of the minus-end and thus filaments are still growing. Below the C_c actin does not polymerize. At the C_c actin monomers and filaments are in steady state, and F-actin is preferentially ADP-actin.

A typical phenomenon of the steady state is treadmilling, where monomers add to the plus-end while others dissociate from the minus-end. For actin alone, the equilibrium is a dynamic exchange of monomers between the G-actin and the F-actin pool. This dynamic is modulated by ABPs (actin-binding proteins).

At physiological salt concentration in the presence of ATP, two factors should be noted influencing the actin polymerization in the absence of ABPs: the actin-concentration and the state of the bound nucleotide.

Means for the choice of experimental conditions

✓ **buffer control**

Add salt and ATP. Polymerization buffers should contain about 0.1M KCl, an excess of 1mM ATP, and 2mM MgCl_2 . ATP-buffers must be kept on ice should be used within 5 days.

✓ **time control**

Wait 15-30min until polymerization reaches steady state.

✓ **temperature control**

Work at room temperature. Polymerization kinetics is usually monitored at 25°C .

Protocol Section

Protocol 1. Preparation of G-Actin

For the preparation of G- or F-actin, the G-actin powder has to be reconstituted to make a stock solution.

1. To obtain a 1.1mg/ml (26.2 μ M) actin stock solution, add 500 μ l H₂O to

Protocol Section

Protocol 2. Preparation of F-Actin

A frequent use of actin is to prepare filaments for cosedimentation assays.

1. Prepare a G-actin stock solution according to Protocol 1.
2. Add 1ml of H₂O to the tube containing PolyMix, and vortex for ~30secs to fully dissolve the PolyMix (10x PolyMix stock solution).

Protocol Section

Protocol 3. F-Actin Cosedimentation Assays at High Speed

Pelleting of F-actin binding ligands is achieved by spinning the samples at 100.000xg for 1h. Cosedimentation assays are best performed at a final actin concentration of 0.5-1.0mg/ml (11.9-23.8 μ M).

Protocol Section

7. Use 40µl of the each tube to prepare a total for SDS-PAGE sample.
8. Use a sample volume of 200µl for high speed spinning. Adjust the volume in the ultracentrifuge tubes using a balance. Slowly overlay sucrose cushions with the sample.
9. Spin the samples at 100.000xg, 1h, 4°C.
10. Unload the rotor when the run is finished, and process the samples immediately.

Protocol Section

Protocol 4. F-Actin Bundling Assays at Low Speed

Bundling of actin filaments is induced by various cellular proteins. This property of proteins has often been discovered by incubation of recombinant fragments with F-actin, as these sites can be cryptic in the full length protein.

For this assay a simple positive control is very useful. In the presence of 50mM $MgCl_2$ to the F-actin stock, actin bundles (so called paracrystals) are formed.

Protocol Section

5. Take 40µl of the each sample to prepare a total for SDS-PAGE.
6. Use a sample volume of 200µl for low speed spinning.
7. Spin the samples at 15.000xg for 30min at 4°C.
8. Unload the rotor when the run is finished and immediately process the samples.

Troubleshooting Guide

Protein-
concentration of
ligand is low

Low ligand concentrations often require the addition of large volumes to the F-actin stock solution. An easy way to compensate low ligand concentrations is to prepare a G-actin stock of 2 or 3mg/ml, by adding 0.45ml or 0.3ml of H₂O to 1mg of G-actin according to Protocol 1.

Ligand Buffer is
not compatible

Phosphate buffers are not compatible with G-actin and should not be used for F-actin either. The best way to overcome

vortexed etc. while polymerizing. Leave the F-actin solution for another 30min at RT and perform the assay.

Bands below actin occur in SDS-

A visible proteolytic actin breakdown product of ~40kD may occur, if actin has been at room temperature for a longer

Toolkit Ordering Information

Product	Description	Cat. #	Size
Actin-Toolkit F-Actin Binding (α -skeletal muscle actin or α -cardiac actin)	Determination and quantification of F-actin binding or bundling by ligands in solution.	8010-01	4x1.0mg Rabbit skeletal muscle actin
Actin-Toolkit G-Actin Binding (α -skeletal muscle actin or α -cardiac actin)	Determination and quantification of ligands binding to ActinBeads in solution.	8020-01	4x250 μ l Rabbit skeletal muscle actin
Actin-Toolkit Fluorometry (α -skeletal muscle actin)	Kinetic measurements of actin dynamics in solution based on pyrenyl fluorescence of actin.	8030-01	8x1.0mg Rabbit skeletal muscle actin
Actin-Toolkit TIRFM (α -skeletal muscle actin or α -cardiac actin)	Single molecule imaging of ligands interacting with ATTO-fluorescent G- or F-actin.	8093-01	4x100 μ g Rabbit skeletal muscle actin
Actin-Toolkit Crystallography (α -skeletal muscle actin)	Co-crystallization of a ligand with non-polymerizable, native G-actin.	8050-01	8x1.0mg Rabbit skeletal muscle actin
Actin-Toolkit SPR (α -skeletal muscle actin or α -cardiac actin)	A unique method to analyze ligand interactions with actin filaments by surface plasmon resonance.	8090-01	4x250 μ g Rabbit skeletal muscle actin
Actin-Toolkit ELISA (α -skeletal muscle actin or α -cardiac actin)	Molecular imaging analysis of ligands bound to monomeric actin, filaments or networks by TEM.	8070-01	4x0.5mg Rabbit skeletal muscle actin
Actin-Toolkit Fluorescence Microscopy (α -skeletal muscle actin or α -cardiac actin)	Identification of ligands bundling actin filaments or forming filament networks by using ATTO-fluorescent actin.	8080-01	4x100 μ g Rabbit skeletal muscle actin

For further Toolkits and more product information we welcome your visit on

www.hypermol.com

Hypermol Services (24/7)

CustomerService

mail : cusserv@hypermol.com

fon : ++49 52 19 87 62 28

TechnicalService

mail : techserv@hypermol.com

fon : ++49 52 19 87 62 30

Order

mail : order@hypermol.com

fon : ++49 52 19 87 62 29

fax : ++49 52 19 87 62 27

We welcome your requests and questions!

Actin-Toolkits and other products from Hypermol are exclusively and directly distributed by:

HYPERMOL
Schulenburger Landstr. 132A
D-30165 Hannover
Germany

mail: hypermol@hypermol.com

fon : ++49 52 19 87 62 26

fax : ++49 52 19 87 62 27

Mail Address: Hypermol | P.O. Box | 201025 D-33549 Bielefeld | Germany

Copyright by HYPERMOL, Germany. License is granted to copy for internal use only.
