

DATASHEET

CapZ

Skeletal muscle, rabbit
(~92% purity)

For Research Use Only.
Not for Use in Diagnostic Processes.

Quantity: 2x25µg
Cat. #: 8320-01

Product Description

The subunits of the heterodimeric capping protein CapZ possess a molecular mass of 36kDa (α -subunit) and 32kDa (β -subunit). CapZ binds with high affinity in a Ca^{++} -independent manner to the barbed end of actin filaments.

The protein concentration of CapZ was determined at UV280nm $0.1\% = 1.234$ and the purity by scanning densitometry of Coomassie G-250 stained SDS-Gels. The lyophilized CapZ contains 100mM KCl, 20mM Hepes pH 7.4, 1mM DTT and 3% disaccharides, when reconstituted with 50µl H_2O .

Determination of capping activity: F-actin (Cat. #:8101-01) containing 20% pyrene actin (Cat. #: 8112-01) was incubated with CapZ and diluted below the critical concentration. Depolymerisation was followed by fluorometry and compared to an actin reference without CapZ. Complete capping was obtained at a molar ratio of 1:100 (CapZ:Actin). All proteins were used as supplied, i.e. without dialysis. Actin stock solutions were pre-spun 15min, 15.000xg, 4°C.

Product Handling

Preparation of a working stock

Add 50µl H_2O to the tube with CapZ to obtain a 1mg/ml stock. Vortex mildly and allow the protein to rehydrate for 2min and vortex again. Quickly spin the tube up to ~3.000xg for seconds to collect the solution at the bottom. Depending on the purpose this stock solution is ready-to-use. For critical assays we recommend buffer exchange against 100mM KCl, 20mM Hepes pH 7.4, 1mM DTT.

Product Storage and Stability

Store the product at -70° upon arrival, where it will be stable for at least 6 months. Once dissolved, CapZ is kept on ice and should be used within 5 days. Avoid refreezing. If aliquots shall be prepared for storage, add glycerol to a final concentration of 50% and store at -20°C .

HYPERMOL EK WARRANTS THAT ITS PRODUCTS ARE CONFORM TO THE INFORMATION CONTAINED IN THIS PUBLICATION. THE PURCHASERS MUST DETERMINE THE SUITABILITY OF THE PRODUCT(S) FOR THEIR PARTICULAR USE. ADDITIONAL TERMS AND CONDITIONS MAY APPLY (SEE REVERSE SIDE OF THE INVOICE). HYPERMOL BRANDS ARE SOLD EXCLUSIVELY BY HYPERMOL EK, GERMANY AND AUTHORIZED DISTRIBUTORS.

For product inquiries please contact:

cusserv@hypermol.com
Fon: +49 (0)521 9876228
Fax: +49 (0)521 9876231
www.hypermol.com